dandy | VISION

Instruction Manual

Contents

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

EN/English

Notice

We recommend that you thoroughly familiarize yourself with this guide to make the most effective use of your system.

The content of this document may change over time. The manufacturer and its affiliates are not responsible for any errors in this document or any incidental damages related to its provision, performance, or use. Reproduction of any part of this publication is prohibited without the prior written permission of the manufacturer.

This document is originally written in English.

All trademarks and registered trademarks belong to their respective owners.

The Dandy Vision is intended for professional use only.

Under U.S. federal law, this device can only be sold by a dentist or on a dentist's order.

If any serious incident occurs in relation to the device, the user must report it to the manufacturer and to the competent authority of its Member State in the European Union.

Conventions

The following special notices highlight important information or potential risks to personnel or equipment.

	WARNING: Indicates safety instructions that must be followed precisely to prevent injury to yourself or others.
	CAUTION: Warns of conditions that could cause severe damage or operational issues.
	NOTE: Provides additional information and helpful tips.


Overview

The scanner is designed to acquire 3D models in the followings:

- Upper jaw
- Lower jaw
- Buccal bite registration

Nomenclature

The scanner consists of the following parts.

The scanner needs to be connected to a computer. Please note that the computer is not included.

[1] Handpiece
[2] Scanner Tip
[3] Scan Button ● Press to start scanning. Press again to stop scanning.
[4] Mode Switch Button ● Press this button to switch between different modes.

[5] Status Indicator ● Rapidly flashing: Scanner is booting. ● Slowly flashing: Scanner is ready but not selected as the active scanner. ● Steady on: Scanner is selected as the active scanner. ● Breathing: Scanner is ready to be connected to the scanning software.
[6] Lens Window
[7] USB Cable (1.8 m, shielded)
[8] USB Type-C Connector
[9] USB Type-A Connector

Safety Information

Intended Use

The Dandy Vision Intraoral Scanner, model DI-1 (hereinafter referred to as "the scanner"), is a digital optical scanning device used to obtain digital impressions of hard and soft tissues, such as teeth, gingiva, and mucous membranes, using oral scanning, for oral restoration and orthodontic treatment of malocclusion.

The Dandy Vision could be used for both adults and children in clinical practice.

Clinical Benefits and Performance Characteristics

The Dandy Vision Intraoral Scanner benefits a dental practice by enabling practitioners to acquire digital impressions with the quality and accuracy required for digital CAD/CAM dental applications. The actual performance of the device is dependent on the user's training and operating execution. The user is solely responsible for the accuracy, completeness, and adequacy of the acquired data.

Contraindications

None.

Warnings and Safety Instructions

DANGER OF ELECTRIC SHOCK: This is an electrical unit. **DO NOT** expose it to water spray. Such action can cause an electric shock or a malfunction of the unit.

CAUTION: All known residual risks, contraindications, or undesirable side effects are listed in this guide. If any serious incident occurs in relation to the device, you must report it to the manufacturer and to the competent authority of your Member State in the European Union.

WARNING

Dandy Vision

- Before using the scanner, you must read and understand all safety information.
- Use the scanner only within hospitals and professional healthcare facilities. Do not operate it near high-frequency surgical equipment or inside the RF-shielded room of an ME system for magnetic resonance imaging, where electromagnetic disturbances are intense.
- Before operating the scanner, inspect the scanner and its accessories for any rough surfaces, sharp edges, or protrusions that could pose a safety hazard.

- You are responsible for the scanner's operation and maintenance. Training is required to use the scanner safely.
- Do not place objects within the scanner's operational field.
- Ensure the scanner is turned off when not in use.
- Do not use the scanner in oxygen-rich environments. It is not intended for use with flammable anesthetics or agents.
- Avoid pulling or twisting the cable.
- Do not drop the scanner or its accessories.
- Do not attempt to sterilize the scanner.
- Keep the scanner away from water sprays, and do not immerse it in water or disinfectants.
- Avoid subjecting the scanner to heavy vibrations.
- Avoid exposing the scanner to ultraviolet radiation directly. The scanner is not designed for ultraviolet disinfection.
- Avoid staring directly at the LED emission window.
- Avoid removing the cover from any part of the scanner. The device contains no user-serviceable components. If repairs are needed, please contact a qualified service technician authorized by the manufacturer.
- Do not use any cables other than those provided with the scanner. Using different cables may damage the scanner and compromise its safety features and EMC performance.
- Keep equipment that does not comply with IEC-60601 standards at least 1.5 meters from the patient.
- If the equipment malfunctions, switch it off, post an "Out of Service" notice, and contact your local distributor or authorized service provider.
- The use of components, accessories, cables, or spare parts not specified or supplied by the manufacturer of this equipment can compromise the scanner's safety features. This may lead to higher electromagnetic emissions, reduced electromagnetic immunity, and improper equipment functioning.
- Do not modify this equipment in any way.
- Do not connect additional power strips or extension cords to the system.
- The applied part may reach temperatures up to 43 °C (109.4 °F). To prevent overheating, avoid prolonged use.
- To completely disconnect the device from the main power supply, unplug the USB connector from the USB port.
- DO NOT maintain or service this equipment while it is in use with the patient.
- Connection of the PEMS (Programmable Electrical Medical System) to an IT NETWORK that includes other equipment could result in risks to patients, operators, or third parties. The responsible organization should identify, analyze, evaluate, and control these risks.

- When using the device on patients with oral mucosal diseases, mental illnesses, severe respiratory diseases, asthma, Parkinson's disease, attention deficit hyperactivity disorder (ADHD), and epilepsy, use extreme caution.
- Patients with moderate or severe opening limitation should use it with caution.

Computer

- DO NOT place any equipment which does not comply with IEC 60601-1 in the immediate vicinity of the patient. Leave at least 1.5 meters distance between the patient and the equipment.
- The scanner is only intended to be connected to a computer that is at least IEC 60950 / IEC 62368, or equivalent standards certified. Connecting the scanner to other equipment may be hazardous.
- Refer to your computer's installation guide for details on the data processing system, computer, and monitor. Ensure there is enough space around the computer to allow for proper ventilation.
- Place the screen so that light reflections from both internal and external sources are minimized to achieve the best image quality and ensure comfortable viewing.
- The equipment under test (EUT) is intended for use only with a computer certified to EN 55032 or an equivalent standard. The computer used during type testing was part of the test setup and is outside the scope of this evaluation.

Disposal

This device contains materials and chemicals commonly used in electrical and electronic products. Improper disposal at the end of its life can lead to environmental pollution. Therefore, do not discard it with regular household waste. Please take it to a designated electronic waste disposal or recycling facility. For more information on how to dispose of electronic waste, please get in touch with your local authorities. Dispose of the Scanner Tips according to standard operating procedures or local regulations for the disposal of contaminated medical waste. For additional Scanner Tips, contact your dealer.

Cleaning, Disinfecting, Sterilizing

Regularly perform the maintenance tasks outlined below for your scanner and accessories.

Carefully follow the instructions when preparing the scanner for use to ensure the highest level of hygienic safety for patients.

To ensure maximum hygienic safety for the patient and to minimize the risk of cross-contamination, carefully perform the following maintenance activities on your scanner and accessories.

After each patient:

- Clean and disinfect the scanner. (See "[Clean and Disinfect the Scanner](#)")
- Clean the Scanner Tip, and then perform Autoclave Sterilization (See "[Clean the Scanner Tip](#)" and "[Sterilize the Scanner Tip](#)").

There are 2 product models of Scanner Tip:

Product Model	Size	UDI-DI	Manual Cleaning	134°C Sterilization	132°C Sterilization
TPDI1L	Medium	(01)06973993441942	Yes	Yes	Yes
TPDI1S	Small	(01)06973993441959	Yes	Yes	Yes

Clean and Disinfect the Scanner

General Warnings

WARNING

- When handling the scanner, read and follow all warnings and personal protection guidelines in the sanitizer's Safety Data Sheet (SDS).
- Gloves must be worn when cleaning and disinfecting the scanner.
- After each patient, disinfect the scanner using a recommended medium-level disinfectant wipe with tuberculocidal properties.
- Do not use disinfectants containing phenol or iodine-based substances, as they may damage the scanner's surface coating.
- Do not place the scanner in a pressure steam sterilizer or immerse it in water or disinfectant solutions.
- Excessive liquid can harm the scanner.
- Avoid using cotton balls, cloth pieces, or paper towels soaked in disinfectant to sanitize the scanner.

Clean the Scanner

If visible blood or body fluid is on the scanner, clean it before proceeding with disinfection.

To clean the scanner, follow the below steps:

- Dampen a lint-free cloth with warm water (do not soak).
- Use the dampened lint-free cloth to wipe away blood or body fluids.

Disinfect the Scanner

The scanner must be thoroughly sanitized after use for each patient.

To properly disinfect the scanner, make sure to follow the disinfectant manufacturer's guidelines regarding the required contact time.

CAUTION: The scanner must be cleaned before disinfection if it is visibly soiled. Refer to "[Clean the Scanner](#)."

To disinfect the scanner, proceed with the following steps:

- Detach the Scanner Tip.
- Clear any visible debris (refer to "[Clean the scanner](#)").
- Use an intermediate-level disinfectant wipe that is commercially available. Be sure to follow the manufacturer's instructions for the necessary contact time.

Recommended disinfectant wipes: CaviWipes

WARNING: Using an unapproved disinfectant may damage the scanner.

- Thoroughly wipe down all scanner surfaces, ensuring no liquid enters any gaps, air outlets, or pinholes.

WARNING: Do not rinse the scanner.

- Let it air dry thoroughly.
- Once dry, use a clean, lint-free cloth dampened with water to remove any remaining disinfectant from the scanner's surface.

Clean and Sterilize the Scanner Tip

WARNING

- Always wear gloves when handling a contaminated Scanner Tip.

- Review and follow the safety and personal protection instructions in the detergent manufacturer's Safety Data Sheet (SDS).
- Avoid soaking the Scanner Tip in disinfectant for extended periods.
- Ensure the scanner tip is completely dry before mounting it to the scanner.
- Do not use an ultrasonic cleaner to clean the Scanner Tip.
- Do not soak the Scanner Tip with alcohol-based disinfectants.

Clean the Scanner Tip

To manually Clean the scanner Tip, follow these steps:

- Rinse the scanner tip under running water to remove excess debris (approximately 2 minutes).
- Brush all surfaces with an enzymatic detergent solution, like Metrex EmPower.
- Rinse thoroughly under fresh, running water for a full 2 minutes.
- Inspect the scanner tip. If it remains dirty, go through the steps above again.
- Clean the mirror inside the Scanner Tip using a lint-free cloth or lens tissue to remove dust.
- Dry the scanner tip's mirror carefully with a lens tissue or lint-free cloth.

Sterilize the Scanner Tip

Scanner Tips shipped non-sterile. You must sterilize them before use.

CAUTION: For TPDI1L and TPDI1S, if you limit the exposure time at 134°C to not more than 6 minutes, you can autoclave the Scanner Tip up to 180 cycles.

To sterilize the Scanner Tip, follow these steps:

- Put the Scanner Tip into a sealed steam sterilization pouch.
- Place the Scanner Tip in a steam autoclave for sterilization.
 - ◆ Exposure temperature should be set as 134°C.
 - ◆ Exposure time should exceed 3 minutes.
 - ◆ Exposure time should not exceed 6 minutes.

CAUTION: For TPDI1L and TPDI1S, if you limit the exposure time to 4 minutes at 132 °C, you can autoclave the Scanner Tip up to 180 cycles.

To sterilize the Scanner Tip, follow these steps:

- Put the Scanner Tip into a sealed steam sterilization pouch.
- Place the Scanner Tip in a steam autoclave for sterilization.
 - ◆ Exposure temperature should be set as 132°C.
 - ◆ Exposure time: 4 minutes.
 - ◆ Drying time: Follow the autoclave manufacturer's recommendation.

Precautions before Use

Perform the following activities on your product and accessories before use.

Visually Inspecting the Scanner for Damage

Visually inspect the scanner for damage or signs of deterioration by doing the following:

- Inspect the scanner's Lens Window.
- Inspect around the scanner buttons and contact points.

If damage is noted, do not use the scanner, and contact your representative.

Visually Inspecting the Scanner Tips

Visually inspect the Scanner Tips for signs of deterioration by doing the following:

- Verify that the Scanner Tip is not damaged, and its components are not detached.
- Verify that the Scanner Tip mirror does not have any smudges or scratches on it.
- If deterioration is noted, replace the Scanner Tip.

WARNING

The Lens Window on the scanner is a delicate optical component. Mount the lens window protective cover to protect the Lens Window from damage and dirt when the scanner is not in use.


The mirror in the Scanner Tip is a delicate optical component. Its clean and undamaged surface is critical to scan quality.

If you see poor scan quality in the software, clean the mirror of the Scanner Tip and the Lens Window using a microfiber cleaning swab, applying ethanol that is free of impurities.


Setting Up the Scanner

To set up the scanner, follow these steps:

- Go to the official support website and download the installation file and instructions according to the product model.
- Double-click the software installation file.
- Choose a language from the drop-down list and click Ok to install.
- Follow the instructions on the screen to complete the installation.
- Insert the USB Type-C Connector. Then push up the cable thumb push sleeve to secure the USB Type-C Connector.


- To disconnect the USB Type-C Connector, use thumb and index finger to pinch the area with red circle and pull the cable thumb push sleeve off. Then detach the USB Type-C Connector from its socket.

- Remove the lens window protective cover on the Lens Window.

- Firmly mount the Scanner Tip onto the Lens Window.

- Insert the Type-A Connector to a USB 3.0 Type-A port on the computer.

CAUTION: Make sure the scanner is connected to the USB 3.0 port. If connected to a USB 2.0 port, the scanner may not work properly.

- The scanner will automatically turn on, the status indicator will flash rapidly, then the Status Indicator will stay on once it's connected to the computer and ready to scan.

NOTE: When the scanner is connected to the software, the indicator stays in solid on.

Getting Started

Accessing the Software User Interface

To access the software user interface, follow these steps:

- Double-click the software icon on your desktop.
- Wait for the software to open.
- Follow the on-screen instructions to use the software.

Preparing the Teeth

- If there is a preparation area, retract the gingiva by gingival restriction cords. And extract the cords just before scanning the preparation.
- Before starting the scan, dry the teeth thoroughly.
- During the scan, re-dry the teeth moderately.

Preparing the Scanner


The Scanner Tip attached to the scanner provides sanitary shield for patients. Always disinfect the body of the scanner and perform sterilization on the Scanner Tip after each use.

CAUTION: Scanner Tip received from the manufacturer are NOT sterilized. You must sterilize them before the first use.


CAUTION: Avoid any liquid from leaking into the air outlet near the Lens Window or the air inlet at the rear of the scanner (see the figure below), otherwise the scanner may be damaged.

To prepare the scanner, follow these steps:

- Make sure the Lens Window at the base of the scanner is clean by wiping it with a moist, lint-free cloth or lens tissue.
- Slide the Scanner Tip onto the scanner as shown below.

Starting Scanning

To start scanning, place the Scanner Tip on the surface of the tooth to stabilize the scanner and press the Scan button. Wait until a 3D model appears on the 3D model display screen, and then slowly move it along the arch at 0-


5mm from the teeth.

Scanning Protocol

The recommended scanning protocol consists of 3 sweeps: occlusal, lingual and buccal to ensure good data coverage of all surfaces.

It is recommended to start the first sweep from the occlusal surface, you should start with the first molar.

The second sweep can scan both the lingual and buccal sides, and the third scan covers the opposite side of the second sweep.

Troubleshooting

DI-1 Troubleshooting Instructions

Problem Description	Action
No connection between the scanner and the computer.	Detach the USB Type-C Connector from the Handpiece, and insert again. If the problem is not solved, replace the cable.
The 3D model shows discrepancies and overlapping areas.	Perform the scan using the suggested scanning method. Use the Cut tool to eliminate mismatched data and excess tissue, then conduct a rescan.
A decline in precision or poorly stitched images may be noticed during acquisition.	Clean the lens window at the scanner's base using a moist, lint-free cloth or lens tissue to remove dust or water stains. Ensure the Scanner Tip is securely attached.
The Scanner Tip is installed but not detected. The Scanner Tip not detected icon is displayed at the bottom-right of the interface.	After removing the scanner tip, wipe the Scanner Tip Detection Sensor to ensure it is free from dirt or debris. Reattach the Scanner Tip and verify that it is securely and firmly connected to the scanner.
The inner surface of the scanner's lens window at its base shows signs of fogging.	Attach a fully dry Scanner Tip to the scanner, place the scanner in the desk, and wait for the fogging to clear. If the fogging persists after 24 hours, contact your local service provider for help. Ensure the Scanner Tip is completely dry before attaching it, and avoid using disinfectant-soaked cloths to clean the scanner.
The scanner does not cast light.	Ensure that the scanner is connected to a USB 3.0 port. Ensure that your computer meets the Minimum System Requirements. Close the software, and then launch the software again.
After stopping scanning, the scanner is disconnected and cannot be automatically reconnected.	Try the following steps in order. Check if the scanner is turned off. Press any button on the scanner to turn on the scanner again. The scanner startup process may take several seconds. Disconnect the USB connector of the scanner from the computer, and then reconnect it to the computer. The scanner will automatically power on and connect.

Regulatory Information

The device complies with the following regulations:

MDR: (EU) 2017/745 Medical Device Regulation, Class I.
Medical Devices Regulations (Canada).
RoHS: Restriction of the Use of Certain Hazardous Substances in Electrical and Electronic Equipment, 2011/65/EU Annex II and its amendment Directive (EU) 2015/863
RED: Directive 2014/53/EU The Radio Equipment Directive
FCC: Part 15 of The Federal Communications Commission Rules
ISED: Innovation, Science and Economic Development Canada

Compliance with European and International Standards

EN 60601-1 / IEC 60601-1: Medical Electrical Equipment, Part 1: General Requirements for Basic Safety and Essential Performance

ANSI / AAMI ES 60601-1: Medical Electrical Equipment - Part 1: General requirements for basic safety and essential performance

CAN / CSA-C22.2 No. 60601-1: Medical Electrical Equipment - Part 1: General requirements for basic safety and essential performance

EN 60601-1-2 / IEC 60601-1-2: Medical Electrical Equipment, Part 1-2: General requirements for basic safety and essential performance – Collateral Standard: Electromagnetic disturbances – Requirements and tests

EN 80601-2-60 / IEC 80601-2-60: Medical electrical equipment — Part 2-60: Particular requirements for the basic safety and essential performance of dental equipment

EN 62471 / IEC 62471: Photobiological safety of lamps and lamp systems: Equipment classification, requirements, and User's Guide

EN ISO 17664: Processing of health care products – Information to be provided by the medical device manufacturer for the processing of medical devices

EN ISO 17665-1: Sterilization of health care products — Moist heat — Part 1: Requirements for the development, validation, and routine control of a sterilization process for medical devices

EN 60601-1-6 / IEC 60601-1-6: Medical Electrical Equipment, Part 1-6: General requirements for basic safety and essential performance - Collateral Standard: Usability

EN 62366-1/ IEC 62366-1: Medical devices - Application of usability engineering to medical devices

EN 62304 / IEC 62304: Medical device software - Software life cycle Processes

EN ISO 10993: Biological evaluation of medical devices

EN ISO 14971: Medical devices - Application of risk management to medical devices

EN ISO 15223-1: Medical devices - Symbols to be used with medical device labels, labeling and information to be supplied - Part 1: General requirements

EN ISO 20417: Medical devices — Information to be supplied by the manufacturer

ISO 9687: Dentistry - Graphical symbols for dental equipment

AAMI TIR 12: Designing, testing, and labeling reusable medical devices for reprocessing in health care facilities: A guide for medical device manufacturers

AAMI TIR 30: A compendium of processes, materials, test methods, and acceptance criteria for cleaning reusable medical devices

Classification in Accordance with EN/ IEC 60601-1

Type of protection against electric shock: N/A

Degree of protection against electric shock: Type BF Applied Part. Scanner tip is considered as applied part.

Mode of operation: Continuous operation

Flammable anesthetics: Not suitable for use in the presence of flammable anesthetics or a mixture of flammable anesthetics with air or oxygen or nitrous oxide.

Conformity with EN/ IEC 60601-1-2

IEC 60601-1-2 EMC requirements and tests, Medical Electrical Equipment including CISPR 11 Group 1, Class B.

Electromagnetic Compatibility Precautions

Medical electrical devices need specific measures to ensure electromagnetic compatibility (EMC).

Medical devices should be set up and used following the EMC guidelines detailed in this document.

Even equipment that meets CISPR emission standards may disrupt communication with the DI-1.

WARNING: To maintain the DI-1's optimal performance, ensure that all portable and mobile RF communication devices—including peripherals like antenna cables and external antennas—are used at a distance of no less than 30 cm (12 inches) from any part of the system, including the cables specified by the manufacturer. Using such equipment closer than this recommended distance may result in decreased performance of the DI-1.

Guidance and Manufacturer's Declarations

Guidance and Manufacturer's Declaration - Electromagnetic Emissions

The device is intended for use in the electromagnetic environment specified below. The customer or user of the DI-1 should assure that it is used in such an environment.

Emissions Test	Compliance	Electromagnetic Environment – Guidance
RF Emissions CISPR 11	Group 1 Class B	The DI-1 uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.

Guidance and Manufacturer's Declaration – Electromagnetic Immunity for Equipment and Systems

The device should be operated only within the electromagnetic environments described below for optimal performance. Customers or users must ensure that the device is used under these specified conditions.

Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment – Guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±8 kV contact ±15 kV air	±8 kV contact ±15 kV air	For optimal conditions, floors should be constructed from wood, concrete, or ceramic tile. If synthetic materials are used to cover the floors, please ensure that the relative humidity is maintained at a minimum of 30%.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	30 A/m	30 A/m	For optimal conditions, power-frequency magnetic fields should be maintained at levels commonly found in standard commercial or hospital environments.
Radiated RF IEC 61000-4-3	3 V/m 80MHz – 2.7GHz	3 V/m 80MHz – 2.7GHz	In a professional healthcare setting: Warning: To prevent potential interference, portable RF communications devices (including accessories like antenna cables and external antennas) must be kept at least 30 cm (12 inches) away from any part of the DI-1, including any cables specified by the manufacturer. Failure to maintain this distance could reduce the equipment's performance.

The device meets the immunity requirements for proximity fields from RF wireless communication devices, in line with the test levels outlined below according to the IEC 60601-1-2 standard. It is the customer or user's responsibility to ensure that the device is operated in a suitable environment.

Test Frequency (MHz)	Band (MHz)	Immunity Test Levels
385	380 – 390	Pulse modulation 18Hz, 27V/m
450	430 – 470	FM, ± 5 kHz deviation, 1 kHz sine, 28V/m
710	704 – 787	Pulse modulation 217Hz, 9V/m
745		
780		
810	800 – 960	Pulse modulation 18Hz, 28V/m
870		
930		
1720	1700 – 1990	Pulse modulation 217Hz, 28V/m
1845		
1970		
2450	2400 – 2570	Pulse modulation 217Hz, 28V/m
5240	5100 – 5800	Pulse modulation 217Hz, 9V/m
5500		
5785		

Accessories

Using cables or accessories other than those specified—unless they are manufacturer-approved replacements for internal components—could lead to increased emissions or reduced immunity of the medical device.

Other Equipment

WARNING: Avoid placing this equipment next to or stacking it with other devices, as this may cause improper functioning. If such placement is necessary, monitor this equipment and the other devices to ensure they operate correctly.

Technical Specification

Product Model

DI-1

DI-1 Technical Specifications

Components	Technical Specifications
Weight	Handpiece: 175g Tip Medium - Sterilization (TPDI1L): 14g Tip Small - Sterilization (TPDI1S): 12g
Color	3D Full Color
Connectivity	USB 3.0
Power Source	USB 3.0 5V, 900mA
Field of View	Tip Medium - Sterilization (TPDI1L): 16mm x 14mm Tip Small - Sterilization (TPDI1S): 12mm x 12mm
Minimum System Requirements of Laptop	Processor: 11th Generation Intel® Core™ i5-11400H or AMD Ryzen™ 7 5700U Memory: 16 GB Graphics Card: NVIDIA GeForce GTX 2060(Notebook) 6GB Disk: 512 GB SSD Operating System: Win 10 (build 18362+) / Win 11, 64 bit Others: USB 3.0 Port Optional: Touch Screen
Minimum System Requirements of Desktop	Processor: 10th Generation Intel® Core™ i5-10600 or AMD Ryzen™ 5 3600 Memory: 16 GB Graphics Card: NVIDIA GeForce GTX 1660 Ti 6GB Disk: 512 GB SSD Operating System: Win 10 (build 18362+) / Win 11, 64 bit Others: USB 3.0 Port Optional: Touch Screen

CAUTION: It is MANDATORY to check that your system configuration is compatible with the computer system requirements for the DI-1 software.

Environmental Requirements

Components	Environmental Requirements
Operating Temperature	15°C ~ 30°C
Transport and Storage Temperature	-10°C ~ 50°C
Operating Relative Humidity	10% ~ 65% RH
Transportation and Storage Relative Humidity	10% ~ 95% RH
Operating Atmospheric Pressure	70 ~ 106 KPa
Transportation and Storage Atmospheric Pressure	60 ~ 106 KPa

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AVERTISSEMENT

Dandy Vision

- Before using the scanner, you must read and understand all safety information.
- Use the scanner only within hospitals and professional healthcare facilities. Do not operate it near high-frequency surgical equipment or inside the RF-shielded room of an ME system for magnetic resonance imaging, where electromagnetic disturbances are intense.

- Before operating the scanner, inspect the scanner and its accessories for any rough surfaces, sharp edges, or protrusions that could pose a safety hazard.
- You are responsible for the scanner's operation and maintenance. Training is required to use the scanner safely.
- Do not place objects within the scanner's operational field.
- Ensure the scanner is turned off when not in use.
- Do not use the scanner in oxygen-rich environments. It is not intended for use with flammable anesthetics or agents.
- Avoid pulling or twisting the cable.
- Do not drop the scanner or its accessories.
- Do not attempt to sterilize the scanner.
- Keep the scanner away from water sprays, and do not immerse it in water or disinfectants.
- Avoid subjecting the scanner to heavy vibrations.
- Avoid exposing the scanner to ultraviolet radiation directly. The scanner is not designed for ultraviolet disinfection.
- Avoid staring directly at the LED emission window.
- Avoid removing the cover from any part of the scanner. The device contains no user-serviceable components. If repairs are needed, please contact a qualified service technician authorized by the manufacturer.
- Do not use any cables other than those provided with the scanner. Using different cables may damage the scanner and compromise its safety features and EMC performance.
- Keep equipment that does not comply with IEC-60601 standards at least 1.5 meters from the patient.
- If the equipment malfunctions, switch it off, post an "Out of Service" notice, and contact your local distributor or authorized service provider.
- The use of components, accessories, cables, or spare parts not specified or supplied by the manufacturer of this equipment can compromise the scanner's safety features. This may lead to higher electromagnetic emissions, reduced electromagnetic immunity, and improper equipment functioning.
- Do not modify this equipment in any way.
- Do not connect additional power strips or extension cords to the system.
- The applied part may reach temperatures up to 43 °C (109.4 °F). To prevent overheating, avoid prolonged use.
- To completely disconnect the device from the main power supply, unplug the USB connector from the USB port.
- DO NOT maintain or service this equipment while it is in use with the patient.
- Connection of the PEMS (Programmable Electrical Medical System) to an IT NETWORK that includes other

equipment could result in risks to patients, operators, or third parties. The responsible organization should identify, analyze, evaluate, and control these risks.

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AVERTISSEMENT

- When handling the scanner, read and follow all warnings and personal protection guidelines in the sanitizer's Safety Data Sheet (SDS).
- Gloves must be worn when cleaning and disinfecting the scanner.
- After each patient, disinfect the scanner using a recommended medium-level disinfectant wipe with tuberculocidal properties.

- Do not use disinfectants containing phenol or iodine-based substances, as they may damage the scanner's surface coating.
- Do not place the scanner in a pressure steam sterilizer or immerse it in water or disinfectant solutions.
- Excessive liquid can harm the scanner.
- Avoid using cotton balls, cloth pieces, or paper towels soaked in disinfectant to sanitize the scanner.

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AVERTISSEMENT : Using an unapproved disinfectant may damage the scanner.

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AVERTISSEMENT

- Always wear gloves when handling a contaminated Scanner Tip.
- Review and follow the safety and personal protection instructions in the detergent manufacturer's Safety Data Sheet (SDS).
- Avoid soaking the Scanner Tip in disinfectant for extended periods.
- Ensure the scanner tip is completely dry before mounting it to the scanner.
- Do not use an ultrasonic cleaner to clean the Scanner Tip.
- Do not soak the Scanner Tip with alcohol-based disinfectants.

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All rights reserved.

Copyright © 2025 Alliedstar.

The content of this instruction is subject to change.

Alliedstar and its subsidiaries shall not be liable for any indirect damages arising from the preparation, presentation, or use of the information contained in this instruction, including any errors. No part of this instruction can be reproduced without the permission of Alliedstar.

Trademarks

All trademarks and registered trademarks are the property of their respective owners.

Lifetime: 6 years

Manufactured for Dandy

1320 N 300 W

Lehi, UT 84043

Alliedstar Medical Equipment Co., Ltd.

No. 1-4, Floor 3, Unit 2, Building 3, No.222, West third section, Waihuan Road, Yanjiang District, Ziyang, 641300, Sichuan, P.R. China

Etkon GmbH

Lochhamer Schlag 6, 82166 Gräfelfing, Germany

www.allied-star.com

Part number: DD0595

Released Date: 2025/12

Rev: 1.1